

DEFB125-defensin, beta 125 Gene

Reference: PTXBC128195

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DEFB-25; beta-defensin 125; beta defensin 25; defensin beta 125		

For research use only.

GENE DESCRIPTION	defensin, beta 125
GENE ID	245938
ORIGIN SPECIES	Human
ACCESSIONS	BC128195
VECTOR	pDONR223
NCBI SYMBOL	DEFB125
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatctctgatctgacctcattatctgtgggtgctaactcgggtgaccaaaggtagctttgaaccccaaaatgttgaagaataatgtaggacattgcagaagacgatgtttagatactgaaagggtacatactcttttaggaacaagctatcatgctgcatttctataatcacatgaat
atactcgacgaccagcatttctgtgtatcacctagaggataaacattggattatagtgtgtggactcttttactggttcccgatctatgttgaatgatctgataacatttgacacaactaaatttgagaaacatgacacctgagaccaatactctgagactactatgccaccatccgagg
ccactactcccgagactactatgccaccatctgagactgctacttccgagactatgccaccacttctcagacagctcttactcataataa