

C21orf129-chromosome 21 open reading frame 129 Gene

Reference: PTXBC128198

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C21orf129; long intergenic non-protein coding RNA 479		
For research use only.			

GENE DESCRIPTION	chromosome 21 open reading frame 129
GENE ID	150135
ORIGIN SPECIES	Human
ACCESSIONS	BC128198
VECTOR	pDONR223
NCBI SYMBOL	C21orf129
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacggaggcagcctccgagccagcccagcagccatggacggaggagccttggagccagccagcagctctctcactggaggggtggacgggacaggaccgctgtcattccagctggaggaggccaggagtctcagctggatgcagtcacctgacctggagagttca
aagtgccataacacaggaaggctcctggaggggatggcctgggatgggtaggacccgggttccacctcaggaggcctctggccatggagtgtaaggagcgtgctaccacaagcctggatcctggcgtgctccgggctctcctcaagccctggcctcatcagacctcaaggct
ccctagacaggggttaattattcacgtggcctgggacccaggggttcgtttgtttgtcttccgacacagtaggggtccacagggaagtcaggatccaccggcaatatttcagtggttatgttgatgtggcaactggcgtgctgaggtatagagaagtaa