

PRSS1-protease, serine, 1 (trypsin 1) Gene

Reference: PTXBC128226

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TRP1; TRY1; TRY4; TRYP1; trypsin-1; TCR V beta 4.1; beta-trypsin; cationic trypsinogen; digestive zymogen; nonfunctional trypsin 1; trypsinogen 1; trypsinogen A; protease, serine 1		

For research use only.

GENE DESCRIPTION	protease, serine, 1 (trypsin 1)
GENE ID	5644
ORIGIN SPECIES	Human
ACCESSIONS	BC128226
VECTOR	pDONR223
NCBI SYMBOL	PRSS1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatccattcctgatccttacccttggcagctgctcttgcgtccccccttgatgatgacaagatcgtgggggctacaactgtgaggagaattctgtccctaccaggtgtccctgaattctggctaccacttctgtgggtgctccctcatcaacgaacagtgggtgtatcagcaggcc
actgctacaagtcccgcatccaggtgagactgggagagcacaacatcgaagtcctggagggaatgagcagttcatcaatgcagccaagatcctccacccccaatagacaggaagactctgaacaatgacatcatgttaatcaagctctctcagtgtagtaataacgcccgcg
tgtccaccatctcttgcaccgccccctccagccaactggcacgaagtgctcatctctggtggggaacactgcgagctctggcgccgactaccagacagctgagtgctgctgagccagctaaagtgtgaagcctctaccctggaaagattaccagcaacat
gttctgtgtgggcttcttggaggaggcaaggattcatgtcaggggtgattctggtggccctgtggtctgcaatggacagctccaaggagttgtctcctgggtgatggctgtgccagaagaacaagcctggagctacaccaaggtctacaactatgtgaatggattaagaacaccatagc
tgccaatagctaa