

CSN1S1-casein alpha s1 Gene

Reference: PTXBC128227

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CASA; CSN1; alpha-S1-casein; casein, alpha; casein alpha s1		

For research use only.

GENE DESCRIPTION	casein alpha s1
GENE ID	1446
ORIGIN SPECIES	Human
ACCESSIONS	BC128227
VECTOR	pDONR223
NCBI SYMBOL	CSN1S1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggcttcattcacctgtctgtggtgtgctcttgcaggcctaaacttccttagataccagaagccttcagaatccatcagagcagtgagcctataccattagaatcaagaggaatacatgaatggatgaacaggcagagaacattctgagagaaaacagactga
tgaatcaaggatactaggaatgagctactcagaactgtgtgtgagcagagcctgagaagatggaatccagcatcagttcatcagtgaggaaatgtctctcagtaagtgtgcggaacagttttgtagactgaacgaatacaaccaacttcagctgcaagctgcccagcagcaaat
tcgcagaatgaatgaaaacagccatgtccagtgcccttccagcagctcaaccaacttgcctaccctatgctgttggactatccacaatcatgcagtatgttccttcccacggtttccgacatctccatccactgctcatgaaaattgaaaaaataacgtcatgctacagtggtg
a