

C21orf67-chromosome 21 open reading frame 67 Gene

Reference: PTXBC128244

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C21orf69; PRED54; chromosome 21 open reading frame 67		
For research use only.			

GENE DESCRIPTION	chromosome 21 open reading frame 67
GENE ID	84536
ORIGIN SPECIES	Human
ACCESSIONS	BC128244
VECTOR	pDONR223
NCBI SYMBOL	C21orf67
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcgaatcggctcattgcagcctcgacgtccaggactcgagcgtatcccgctcagcctctcgagaaggagtgccagaaggctcctctccaggcctcacagagggaacactggctctagttagagacgtcgctcctgccatcaggaccgaagcggcagcgcgaggat
cctgcaggccaagagacaacagcgattacaaccccagtcaccgcttgcagctgacctgcaggtgacgcccctgctgggtgtctcgggtgcagctgcacatcaggggcctctgtggacagaagctctgagtcacacttggccccaaggcgtggagctggagcactgccacgaga
gggatgctgctggctgcggcagcttgcagcccttccctgcaccgagatgtccagctgagcggctgggtgtcacaagttcccggttgggccatcagaggagatcaaaagatcaacccctcccatgggctcctgcctgcctgcaggtgggtttcctgctgtctccagccccaagag
ctccacggactctgccagtagctgtttaaggaggaggaggttttctgatcctctggacatccctgggctggagctatgggtga