



Reference: PTXBC128263

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PAPD2; RBM21; STARPAP; TUTase; URLC6; speckle targeted PIP5K1A-regulated poly(A) polymerase; PAP-associated domain-containing 2; RNA uridylyltransferase; RNA-binding motif protein 21; RNA-binding protein 21; TUTase 6; U6 snRNA-specific terminal uridylyltr		

For research use only.

GENE DESCRIPTION	terminal uridylyl transferase 1, U6 snRNA-specific
GENE ID	64852
ORIGIN SPECIES	Human
ACCESSIONS	BC128263
VECTOR	pDONR223
NCBI SYMBOL	TUT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]