

ADAMTSL3-ADAMTS-like 3 Gene

Reference: PTXBC128389

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ADAMTSL-3; ADAMTS-like protein 3; a disintegrin-like and metalloprotease domain with thrombospondin type I motifs-like 3; punctin-2; ADAMTS like 3		

For research use only.

GENE DESCRIPTION	ADAMTS-like 3
GENE ID	57188
ORIGIN SPECIES	Human
ACCESSIONS	BC128389
VECTOR	pDONR223
NCBI SYMBOL	ADAMTSL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr