



Reference: PTXBC128557

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PIG5; TMEM16J; TP53I5; anoctamin-9; p53-induced gene 5 protein; transmembrane protein 16J; tumor protein p53-inducible protein 5; anoctamin 9		

For research use only.

GENE DESCRIPTION	anoctamin 9
GENE ID	338440
ORIGIN SPECIES	Human
ACCESSIONS	BC128557
VECTOR	pDONR223
NCBI SYMBOL	ANO9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]