

FCGR3B-Fc fragment of IgG, low affinity IIIb, receptor (CD16b) Gene

Reference: PTXBC128562

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD16; CD16b; FCG3; FCGR3; FCR-10; FCRIII; FCRIIIb; low affinity immunoglobulin gamma Fc region receptor III-B; Fc fragment of IgG, low affinity IIIb, receptor (CD16b); Fc gamma receptor IIIb; Fc-gamma receptor IIIb (CD 16); fc-gamma RIII-beta; fc-gamma RI		

For research use only.

GENE DESCRIPTION	Fc fragment of IgG, low affinity IIIb, receptor (CD16b)
GENE ID	2215
ORIGIN SPECIES	Human
ACCESSIONS	BC128562
VECTOR	pDONR223
NCBI SYMBOL	FCGR3B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggcagctgctcctccaactgctctgctactttagttcagctggcatgcggactgaagatctccaaaggctgtggtgttctctggagcctcaatggtacagcgtgcttgagaaggacagtgtagctctgaagtgccaggagcctactccctgaggacaattccacacagtggtttc
acaafgagagcctcatctcaagccaggcctcgagctacttattgacgtgccacagtcacgacagtgtagagtagacagtgccagacaaacctctccacctcagtgacccgggtgacagctagaagtcctatcggctggctgttctccaggccctcgggtgttcaaggagggaag
accctattcactgaggtgtcacagctggaagaacatgctctgcataagggtcacatatttacagaatggcgaagacaggaagtattttcatcataattctgacttccacattccaaaagccacactcaaatagcggctcctacttctgcagggggctgttgggagtaaaatgtctctcag
agactgtgaacatcaccatcactcaaggtttgacagtgtaacacatctcatcttctccacctgggtaccaagtccttctctgtggtgatggtactcttttgcagtggaacacagactatattctctggaagacaacattga