



Reference: PTXBC128575

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PMCA2; PMCA2a; PMCA2i; plasma membrane calcium-transporting ATPase 2; ATPase, Ca++ transporting, plasma membrane 2; plasma membrane Ca(2+)-ATPase; plasma membrane Ca2+ pump 2; plasma membrane calcium ATPase; plasma membrane calcium pump; ATPase plasma mem		

For research use only.

GENE DESCRIPTION	ATPase, Ca++ transporting, plasma membrane 2
GENE ID	491
ORIGIN SPECIES	Human
ACCESSIONS	BC128575
VECTOR	pDONR223
NCBI SYMBOL	ATP2B2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]