



Reference: PTXBC128585

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HUCEP-4; SCAND4; dJ874C20.4; piggyBac transposable element-derived protein 1; cerebral protein 4; piggyBac transposable element derived 1		

For research use only.

GENE DESCRIPTION	piggyBac transposable element derived 1
GENE ID	84547
ORIGIN SPECIES	Human
ACCESSIONS	BC128585
VECTOR	pDONR223
NCBI SYMBOL	PGBD1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]