



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TGIF2LY-TGFB-induced factor homeobox 2-like, Y-linked Gene

Reference: PTXBC128604

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	homeobox protein TGIF2LY; TGIFLY; TGF-beta-induced transcription factor 2-like protein; TGFB-induced factor 2-like protein, Y-linked; TGFB-induced factor 2-like, Y-linked; TGIF-like on the Y; TGFB induced factor homeobox 2 like, Y-linked		

For research use only.

GENE DESCRIPTION	TGFB-induced factor homeobox 2-like, Y-linked
GENE ID	90655
ORIGIN SPECIES	Human
ACCESSIONS	BC128604
VECTOR	pDONR223
NCBI SYMBOL	TGIF2LY
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagagccgcgtgcagacggccggctgagaccacaagccgggtgagaaagacagaccggcgagaccacaagcccgcccaagaacacctcaatcatgtctgagaataacgcagatacaggcagagtcttctccttaccagagacacaagaagaagcgcaggggaaacttgccagccgagtcggttaagatctctccgcgactgagtataagcatcgtgttaaggcctacccttcagaagaagagagaagaacaaatctgtctgcagagaagaccataattgtcttgttgcggattcttaacgtgtttatcaatgtctctgcagacgcattctcccgatagtctcaacacgctgagaacagaccatcatctggccacaaacgcgacaagaatgcccatgcccccactgcagagcaccgaggcgctgtgccggcgaatgcaggccgattgtccagacaatgtacaagcgctgccctgtgcccctgtccaaaggccgcatgtcaagagagaagcaaccagatccggcgatcgcccttcgccaagctaccgggaatg

More info: Please contact us for any question: contact@proteogenix.fr