

SSX1-synovial sarcoma, X breakpoint 1 Gene

Reference: PTXBC128611

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein SSX1; CT5.1; SSRC; cancer/testis antigen 5.1; cancer/testis antigen family 5, member 1; sarcoma, synovial, X-chromosome-related 1; synovial sarcoma, X breakpoint 1; SSX family member 1		

For research use only.

GENE DESCRIPTION	synovial sarcoma, X breakpoint 1
GENE ID	6756
ORIGIN SPECIES	Human
ACCESSIONS	BC128611
VECTOR	pDONR223
NCBI SYMBOL	SSX1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacggagacgacaccttgcacagagaccaggatgatgctaaagcatcagagaagagaagcaaggccttggatgattgacacatacttcttaagaaagagtggaagaaatcagctatgtgtatatgaagagaaactataaggccatgactaaac
taggtttcaaaagtcaccctcccacctttcatgtgtaataaacagccacagactccagggaatgatttgataatgaccataaccgaggatcaggttgaaacatcctcagatgactttcgagggctccacagaatcctccgaagatcatgcccaagaagccagcagaggaacaaatg
attcgaaggaggtgcagaagcatctgcccacaaaacgatgggaacaactgcacccccaggaaaaagcaaatattctgagaagattaagagatctggacccaaaggggaaacatgcctggaccacagactgcgtgagagaaagcagctggtgattatgaagagatcag
cgaccctgaggaagatgacgagtaa