

C10orf125-chromosome 10 open reading frame 125 Gene

Reference: PTXBC129819

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C10orf125; FUCU; FucM; fucose mutarotase; protein fucU homolog		
For research use only.			

GENE DESCRIPTION	chromosome 10 open reading frame 125
GENE ID	282969
ORIGIN SPECIES	Human
ACCESSIONS	BC129819
VECTOR	pDONR223
NCBI SYMBOL	C10orf125
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggcgctgaagggtgtcccccgaactgctgtcccccagctgctctacggcgtggcgggatggggcagggggacgagatggttcttgcggacttgaacttcccgcctccatctgccagtgtgggcccatggagatccgtgcagacggcctgggcatccgcagctcctgga
ggccgtgctgaagctgctgtcccccggacacctatgtggagagtcgggtgcagtcagtgagctgtgtcccgagcacaaggagaggggcctgcagaccccatgtgtggacggagttacgagtcctctacgcagggccgctgtgtgagagccctggcaagatagagaggtttgagtt
ttatgaacgggctaagaaggctttgtgtgttggaacgggatgctga