



Reference: PTXBC129836

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD307b; FCRH2; IFGP4; IRTA4; SPAP1; SPAP1A; SPAP1B; SPAP1C; Fc receptor-like protein 2; IFGP family protein 4; SH2 domain containing phosphatase anchor protein 1; fc receptor homolog 2; immune receptor translocation-associated protein 4; immunoglobulin re		

For research use only.

GENE DESCRIPTION	Fc receptor-like 2
GENE ID	79368
ORIGIN SPECIES	Human
ACCESSIONS	BC129836
VECTOR	pDONR223
NCBI SYMBOL	FCRL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

atgtggaaatgagaaatatacgaacagccatggggccagcgccttctgagcctgcaggttgggagtttgtgaactattgagcatcacaacattttctgatagccctctgtctgtcagttccagttctccctctcctcaacctcaggtctcctggggccagcgtcagttgggga
cctgtcgtgagcttactgtgagggccctgagaggtctcccccaattctgaccaatttatcatgaggatgtacccttgggaacagctcggccccccttgagaggaggggcctcctcaacctctcttctgactcagaaacattctggaaactactctgtgaggccaaacgacctgggggc
ccaagtgcagtgaggagcagtgccactgtccatctcaggaactctatggctatagaaggaacactcatgacagctggagttctctgggagctcttgggtctcctgggttctactgtgttctgttctgtgtatgcctgttccagaagatcagaggaaagtctctgcacatgaacccagaggggc
ttccagggccaaatctcaaggttccacttaacagcccaacccagacatggaggagcgtcagccagttgtatgcaatgctctatagattggatgtggtttattctcaggtctggagcatgcagcagccagaaagctcagcaaacatcagagacactctggagaacaaggactccc
aagtcacatctctctgttgggaagaataaa