



Reference: PTXBC129930

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CABC1; ADCK3; ARCA2; COQ10D4; COQ8; SCAR9; atypical kinase COQ8A, mitochondrial; aarF domain-containing protein kinase 3; atypical kinase ADCK3, mitochondrial; chaperone activity of bc1 complex-like, mitochondrial; chaperone, ABC1 activity of bc1 complex		

GENE DESCRIPTION	chaperone, ABC1 activity of bc1 complex homolog (S. pombe)
GENE ID	56997
ORIGIN SPECIES	Human
ACCESSIONS	BC129930
VECTOR	pDONR223
NCBI SYMBOL	CABC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1