



Antibodies and Proteins

Espace Européen de l'Entreprise
15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

LHFPL3-lipoma HMGIC fusion partner-like 3 Gene

Reference: PTXBC129937

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LHFPL4; lipoma HMGIC fusion partner-like 3 protein; LHFP-like protein 3; lipoma HMGIC fusion partner-like 4; lipoma HMGIC fusion partner-like 3		

For research use only.

GENE DESCRIPTION	lipoma HMGIC fusion partner-like 3
GENE ID	375612
ORIGIN SPECIES	Human
ACCESSIONS	BC129937
VECTOR	pDONR223
NCBI SYMBOL	LHFPL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr