

LOC338588-hypothetical LOC338588 Gene

Reference: PTXBC129979

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC338588
GENE ID	338588
ORIGIN SPECIES	Human
ACCESSIONS	BC129979
VECTOR	pDONR223
NCBI SYMBOL	LOC338588
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctcttgaaaacacggtttgatggagaggtttgtttgttaacctgtatgccgttaatacagaagagccgtctctggcagatatgtgaagaactgccatgaggagcggagaggctacttcgcacacctgcagctcaggtcccgagacgggtttacaggcacacaagccctctccca
agccatctgccaggagcccatgaggaaacacagcacagcgctgcatacagaagccatgggtttgcctatcacatgtctgtttctctctgatccgtggctcacaggctcagttgacccaacctctag