



Reference: PTXBC129993

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT108; otoancorin; cancer/testis antigen 108		
For research use only.			

GENE DESCRIPTION	otoancorin
GENE ID	146183
ORIGIN SPECIES	Human
ACCESSIONS	BC129993
VECTOR	pDONR223
NCBI SYMBOL	OTOA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr