



Reference: PTXBC130009

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CFAP39; CLA2; OPCA; PMCA3; PMCA3a; SCAX1; plasma membrane calcium-transporting ATPase 3; ATPase, Ca++ transporting, plasma membrane 3; cilia and flagella associated protein 39; plasma membrane calcium ATPase; plasma membrane calcium pump; ATPase plasma me		

GENE DESCRIPTION	ATPase, Ca++ transporting, plasma membrane 3
GENE ID	492
ORIGIN SPECIES	Human
ACCESSIONS	BC130009
VECTOR	pDONR223
NCBI SYMBOL	ATP2B3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]