



Reference: PTXBC130294

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	angiomotin; angiomotin p130 isoform; angiomotin p80 isoform		
For research use only.			

GENE DESCRIPTION	angiomotin
GENE ID	154796
ORIGIN SPECIES	Human
ACCESSIONS	BC130294
VECTOR	pDONR223
NCBI SYMBOL	AMOT
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]