

LOC441108-hypothetical gene supported by AK128882 Gene

Reference: PTXBC130299

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	chromosome 5 open reading frame 56		
For research use only.			

GENE DESCRIPTION	hypothetical gene supported by AK128882
GENE ID	441108
ORIGIN SPECIES	Human
ACCESSIONS	BC130299
VECTOR	pDONR223
NCBI SYMBOL	LOC441108
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggactccctggctgctggagagttgaatgccagccaccagccatgggtgccagagttgtagcctattggaggaaaacacaccaaggggaatcttcaaaactcttctccactggctcctgctgtgctcctgtctccatcacctgtataatgcttccctgccatcctgcttccccgtggagtag
agccatgctgtgctgctttacagccaagattctggaagaagaatgtctccttttctcgtcatcaaccactcattattcagcaagtggttttacttgcctgtgctgggcatggcgaatgcagcagtgaaacaaaaccaccagccctgctctcatgcagcttataatcaagtagagagaca
gacataa