

C21orf87-chromosome 21 open reading frame 87 Gene

Reference: PTXBC130301

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	chromosome 21 open reading frame 87
GENE ID	257357
ORIGIN SPECIES	Human
ACCESSIONS	BC130301
VECTOR	pDONR223
NCBI SYMBOL	C21orf87
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcttgatgatccagccgagccgtcagccgagctggcgctgctccccgactcgacccccctggcggaagcgctgcatcttattgtaatttctgccgggctgctgcttgggggatgtggctgggtctcggcagccctcagcggcgagacgtccggcaacgccagcatag
aggccccgccccagggagcggggtcggcggttccccctctgctgggagcgcgcgagccaatcctcactttaaaagcttctctgctcgaagcgaggggtctctcgtgctggggcaagacgcactggccccagacgaccggtcgtgaaatgaccttccgagcgaccgga
gcgggagacctgccccaacccgcgactcggggccctgctggtcccgctagaggacccaggtccacggggggcgcggtga