

SUMO4-SMT3 suppressor of mif two 3 homolog 4 (*S. cerevisiae*) Gene

Reference: PTXBC130305

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IDDM5; SMT3H4; SUMO-4; dJ281H8.4; small ubiquitin-related modifier 4; SMT3 suppressor of mif two 3 homolog 2; SMT3 suppressor of mif two 3 homolog 4; insulin-dependent diabetes mellitus 5; small ubiquitin-like modifier 4 protein; small ubiquitin-like modi		

For research use only.

GENE DESCRIPTION	SMT3 suppressor of mif two 3 homolog 4 (<i>S. cerevisiae</i>)
GENE ID	387082
ORIGIN SPECIES	Human
ACCESSIONS	BC130305
VECTOR	pDONR223
NCBI SYMBOL	SUMO4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaacgaaaagccacagagaagtgcaagctgagaacaacaatcatatttgaaggtggcgggacaggatgggtctgtggtgcagtttaagattaagaggcagacaccacttagtaaaactaatgaaagcctattgtgaaccacggggattgtcagtgaaagcagatcagattcc
gatttgggtgggcaaccaatcagtggaacagacaaacctgcacagttggaaatggaagatgaagatacaattgatgtgtttcaacagcctacgggaggtgtctactga