

TMEM182-transmembrane protein 182 Gene

Reference: PTXBC130306

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 182		
For research use only.			

GENE DESCRIPTION	transmembrane protein 182
GENE ID	130827
ORIGIN SPECIES	Human
ACCESSIONS	BC130306
VECTOR	pDONR223
NCBI SYMBOL	TMEM182
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagactaaatcgtctatcttcttggagctctcttgggtgcttgggggtgttactcttttgggtgcttggatcggattattggctcttgcactgaagtggggagatgttcaggtgaaaagaatatagagaacgtcactttccacatgaagggttcttggaggtgtgtttaaagggtatg
tggagagaatgactccaatatttggagttctggtacaccaatcagccaccgtccaagaactgcacacatgcttacctgtctcgtacccttcatgagaggcagacacactcgacctctatgactctgagttattaccgtgttcttgggcagtcctgatctcctggggtagttgctg
tagtcatcgcaagcttttgcacatctgtgcagcccccttcgcagccatttctctacaaagctgggggaggtcatatattgctgcaggcatctattttcattggtgtgtgctgtatgtcatctgggtccaggcagtgctgacatggaagctaccgaaacatgaaatgaaggactgcct
ggatttcaccccttctgtctgtatggctgtcatcttctcgtgccccagctgggataatttttcttctgtagctggattactatttctgtgtgtggacgcatattcagatacatcactaa