

RNASE9-ribonuclease, RNase A family, 9 (non-active) Gene

Reference: PTXBC130311

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HEL128; RAK1; h461; inactive ribonuclease-like protein 9; ribonuclear enzyme; ribonuclease A K1; ribonuclease, RNase A family, 9 (non-active); ribonuclease A family member 9 (inactive)		

For research use only.

GENE DESCRIPTION	ribonuclease, RNase A family, 9 (non-active)
GENE ID	390443
ORIGIN SPECIES	Human
ACCESSIONS	BC130311
VECTOR	pDONR223
NCBI SYMBOL	RNASE9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgagaactcaccacacacccactgccctctctattgccgagcagctgctgagctggtgcagttcaagaggtggatacagatttgaltcccagaagaagataaaaaagaagaattgaagagtgttggaaaaatttttagtacagggccgccagaccactacaa
agaaaaagtcaaaagacgtctctattgaaactggaatgccactaaatcatatagactgtaacctgaaatcatgggaaaaatgttactacaaacaccgtgggtggcagaacattctctctatgcaatfatgacgagctcaaaaaatctgttacaacagatttgccatgtaag
aatggaattagaaatgtaacaggagcaaaaggctttagaaggagtgattgtaatttaacagaagcatgaaataccagcgtgtaaatacgaatcattataggaagggctacgtcttactctgcatggcaaaatgaaatgcaaaaactattcctcactataaatgatctctgg
agccacctgaacacagaagttctctcagtgaggatggtgtctttgtcataccgccctag