

IL28B-interleukin 28B (interferon, lambda 3) Gene

Reference: PTXBC130314

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IL28B; IFN-lambda-3; IFN-lambda-4; IL-28B; IL-28C; IL28C; interferon lambda-3; cytokine Zcyto22; interferon, lambda 4; interleukin-28B; interleukin-28C; interferon lambda 3		

For research use only.

GENE DESCRIPTION	interleukin 28B (interferon, lambda 3)
GENE ID	282617
ORIGIN SPECIES	Human
ACCESSIONS	BC130314
VECTOR	pDONR223
NCBI SYMBOL	IL28B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccggggactgcatgccagtgcgtgctgatggccgagtgctgacgtgactggagcagttcctgtcgcagggctccggggctctccggatgcaaggggctgccacatagcccagttcaagtcctgtctccacaggagctgcaggccttaagaggccaaagatgcctt
agaagagtcgtctgtgaaggactgcaagtcggtccgctccgctcttcccaggacctgggacctgagggcagctgcaggtgaggagagcggccgcttggaggtgagctggccctgacgtgaaggttctgaggccaccgctgacactgacccagcctgggggatgtcttg
gaccagcccttcacaccctgacccatctctccagctccgggctgtatccagcctcagcccagggcagggccagggacccggggcgccctccaccattggctgcaccggctccaggagggcccaaaaaggagtccttggctgctcagggcctctgacacctcaacctt
ccgctctcacgcgagacctgaattgtgtgcccaggggacctgtgtgtctga