

MT1DP-metallothionein 1D (pseudogene) Gene

Reference: PTXBC130315

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			

For research use only.

GENE DESCRIPTION	metallothionein 1D (pseudogene)
GENE ID	326343
ORIGIN SPECIES	Human
ACCESSIONS	BC130315
VECTOR	pDONR223
NCBI SYMBOL	MT1DP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacctcagctgctcctcgccactgggtgctcctgcacctgtgccagctctgcaaatgcaagagtagcaaatgcacctctgcaagaagaactgctcctcctgccccatgggctgtgccaatgtgccaggctgcacctga