

STH-saitohin Gene

Reference: PTXBC130319

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MAPTIT; saitohin; microtubule-associated protein tau (MAPT) intronic transcript		
For research use only.			

GENE DESCRIPTION	saitohin
GENE ID	246744
ORIGIN SPECIES	Human
ACCESSIONS	BC130319
VECTOR	pDONR223
NCBI SYMBOL	STH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtgagggtggagcccaagtcctcatgcatcttgcagccccacagaactgtgcaggtggccggccctcattgaatgcgggttaatttaactcagcctctgtgtgagtgatgattcaggtgccagagacagaaacccctagcttagcatgggaagtagcttcctgttgaccctgagtt
catctgaggtggcttggaggtgtggcaccatttggcccagttctacagctctgaagagagcagcaggaatgggctgagcaggggaagacaactttccattgaaggccctttcagggccagaactgtccctccaccctgcagctgcctctgcccctgaggggtgagagtc
ggcgacctcatccaagtgtag