



Reference: PTXBC130340

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	glutamine-rich protein 1; glutamine rich 1		
For research use only.			

GENE DESCRIPTION	glutamine-rich 1
GENE ID	54870
ORIGIN SPECIES	Human
ACCESSIONS	BC130340
VECTOR	pDONR223
NCBI SYMBOL	QRICH1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1