



Reference: PTXBC130344

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	gamma-aminobutyric acid receptor subunit rho-1; GABA(A) receptor subunit rho-1; GABA(A) receptor, rho 1; GABA(C) receptor; bA135P14.1 (gamma-aminobutyric acid (GABA) receptor, rho 1); gamma-aminobutyric acid (GABA) A receptor, rho-1; gamma-aminobutyric ac		

GENE DESCRIPTION	gamma-aminobutyric acid (GABA) receptor, rho 1
GENE ID	2569
ORIGIN SPECIES	Human
ACCESSIONS	BC130344
VECTOR	pDONR223
NCBI SYMBOL	GABRR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1