

RNASE12-ribonuclease, RNase A family, 12 (non-active) Gene

Reference: PTXBC130347

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HEL-S-85p; RAI1; epididymis secretory sperm binding protein Li 85p; ribonuclease A I1; ribonuclease, RNase A family, 12 (non-active); ribonuclease-like protein 12; ribonuclease A family member 12 (inactive)		

For research use only.

GENE DESCRIPTION	ribonuclease, RNase A family, 12 (non-active)
GENE ID	493901
ORIGIN SPECIES	Human
ACCESSIONS	BC130347
VECTOR	pDONR223
NCBI SYMBOL	RNASE12
LENGTH	Proteogenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgataataatggtgataatcttctggtgcttctgttctgggaaatgaggtgaatgatgaagcagtgatgcaacttagaaccttgcacgtggactaccctcagaatgacgttccccttctgcaagtgactgcaaccacatgatcataaaagattatcagggaaactgaccacactgtaa
aaaggacatgtcttccatcatgatgagagcctcgaaaaatcaatggtatcttcccaagaaggtgcttgcacaaaccttcggccatttctgcttcagagtgagacaaaagtcaaaatgacagctctgacgctcattgaaaggcacaagataacctgcctgacgtaccactatcccc
acagaggggttctctgacactgtgatgactggagccagatagttccttgctatgtaataa