

BLID-BH3-like motif containing, cell death inducer Gene

Reference: PTXBC130361

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BRCC2; BH3-like motif-containing cell death inducer; breast cancer cell 2; breast cancer cell protein 2; BH3-like motif containing, cell death inducer		

For research use only.

GENE DESCRIPTION	BH3-like motif containing, cell death inducer
GENE ID	414899
ORIGIN SPECIES	Human
ACCESSIONS	BC130361
VECTOR	pDONR223
NCBI SYMBOL	BLID
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgactttgtgctatagaggccaggaatacatcttcttgagatcctagaaatctgagtggtgtctctacacaggatggatagagcgagcctctggcagttccattatccagaggcaaaagcacgcctgccactggaggcgctcttgggtccaacaaagacacatgttgcctaagga
aacagtgctttctctaaagggtacaatcttggtcctctgccatgaagcggaatgttctctggacatgtgcttcagagaccttccatttaaccaggatacaagttacattgttatgcaattcctctgctgaggccctgtaa