

LCE3C-late cornified envelope 3C Gene

Reference: PTXBC130365

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LEP15; SPRL3A; late cornified envelope protein 3C; late envelope protein 15; small proline rich-like (epidermal differentiation complex) 3A; small proline-rich-like epidermal differentiation complex protein 3A; late cornified envelope 3C		

For research use only.

GENE DESCRIPTION	late cornified envelope 3C
GENE ID	353144
ORIGIN SPECIES	Human
ACCESSIONS	BC130365
VECTOR	pDONR223
NCBI SYMBOL	LCE3C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcctgccagcaaaaccagcagcagtcgagcccccctccagttgtccctcaccgaagtgccccaaagagccagcacagtgctgctccaccctcttctgactgtgctctaagctccggggctgtggccccagttctgaaagtggctgctgctgagccaccaggcactca
gggtcccatcaatgccggcgccagagatccaactctgtgacaggggcagtggtcagcaaggcgggggctcctgcccgtggccatggctctgggggctgctgctga