

GAGE12B-G antigen 12B Gene

Reference: PTXBC130374

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GAGE-12B; G antigen 12B/C/D/E; g antigen 12C/D/E; G antigen 12B		
For research use only.			

GENE DESCRIPTION	G antigen 12B
GENE ID	729428
ORIGIN SPECIES	Human
ACCESSIONS	BC130374
VECTOR	pDONR223
NCBI SYMBOL	GAGE12B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttggcgagggaagatcgacattattggcctagaccaaggcgctatgtacagcctctgaaatgattggcctatgcggccgagcagttcagtgatgaagtggaaccagcaacacctgaagaagggaaccagcaactcaatgtcaggatcctgcagctgtcaggaggaggaga
ggatgaggagagcatctgcagggtcaaggccgaagcctgaagctcatagccaggaacagggtcaccacagactgggtgtgagtgatgaagatggctctgatgggcaggagatggacccgccaatccagaggaggtgaaaacgcctgaagaaggtgaaaagcaatcacagtgttaa