



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

KRTAP5-6-keratin associated protein 5-6 Gene

Reference: PTXBC130399

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

KRTAP5.6; keratin-associated protein 5-6; keratin-associated protein 5.6; ultrahigh sulfur keratin-associated protein 5.6; keratin associated protein 5-6

For research use only.

GENE DESCRIPTION	keratin associated protein 5-6
GENE ID	440023
ORIGIN SPECIES	Human
ACCESSIONS	BC130399
VECTOR	pDONR223
NCBI SYMBOL	KRTAP5-6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgctgctgctgctgctggaggctgtggctccggctgtgggggctgtggcttggctgtgggggctgtgggtccagctgctgtgtgccatctgtctgtcgaagccgctgtgctgtgtgtgccacctgttctgcaccagctgtggctctgtggggggctccaaggggtgtgtgtgg
ctctgtggggggctccaaggggtgtgtgtgctgtgctgtgtgggggctccaagggagctgtggctctgtgtgctgtccactgtcagttgtgtcgaagccctgtactgttctcagctgtgtgggtcatctgtcgtccagtcacagctgtgtcgaagccctgtgttccagagccagctgtgtgtcc
catttctcccaagtcacaaactga