

BOLA2-bolA homolog 2 (E. coli) Gene

Reference: PTXBC130401

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BOLA2A; BOLA2B; My016; bolA-like protein 2; BolA-like protein 2 member A; bolA homolog 2; bolA family member 2		

For research use only.

GENE DESCRIPTION	bolA homolog 2 (E. coli)
GENE ID	552900
ORIGIN SPECIES	Human
ACCESSIONS	BC130401
VECTOR	pDONR223
NCBI SYMBOL	BOLA2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaagcgcgaaaagcctggaccgtggaaagcccgctcttgaggcggaagtactgcgttgacgtacgccttagtaaggcggaagtgtttccagcgaagtggctcctgtaaggcagcaagtagcgtggccggcgcccgagctggggtgtgtccctgctgggctgc
cggtccagctggactgcccatggaactcagcgcgaatacctccgcgagaagctgcagcgggacctggaggcggagcaltgtggagggtggaggacacgacctcaaccgtgtcctctgtagcttccgagtcctggtgtgtcggccaagttcgaggggaaccgctgttcagaga
cacaggttctgtacagagtaa