

LOC285679-hypothetical protein LOC285679 Gene

Reference: PTXBC130442

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC285679
GENE ID	285679
ORIGIN SPECIES	Human
ACCESSIONS	BC130442
VECTOR	pDONR223
NCBI SYMBOL	LOC285679
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgtcatcagtcgccgaagacacccgtagagtccttgggtctccatcatcctgagctcctcaagccagagagcctctgtgtccctgaagcacccttcacaccagccacctgcgagcacctatcaccaaaccggaccagctccagaaatccttggggagaccccatgca
ggcaggtggagcaggggtcccttctcagcctggatgtgccaacgctgctggacgcataagacaccaaaccagacactccatggacacgaatatccgacatccacgaacaattgggaatctctga