

LOC158381-ATPase, Class I, type 8B family pseudogene Gene

Reference: PTXBC130448

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	ATPase, Class I, type 8B family pseudogene
GENE ID	158381
ORIGIN SPECIES	Human
ACCESSIONS	BC130448
VECTOR	pDONR223
NCBI SYMBOL	LOC158381
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacagatggcagtgctccagaagccagtcacattagactctacatctctatgcagccaagccctgtatcaaatcagtcactttatcagaatctgtgcatcttctcaattggacagtcacatctgtggacaagcagtcactgaaacagaagatctgcaggctcagtcagacatg
gcacagcagccatctgaagagcaaaagcaagtccttgaataacaaaacaaaaaagaatctctgtttcatgtgcaggaagaaagtgaggacttactaggtttgaatgcggtgtagaaatgtttactgtgtgtacactgttactcagatctacacaattgctcttaaaattacaaagctgatgccg
ctgagaaaatcagaaaagaaatccagtag