

C6orf159-chromosome 6 open reading frame 159 Gene

Reference: PTXBC130460

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf159; SCDO6; dJ237I15.1; protein ripply2; ripply2 homolog; ripply transcriptional repressor 2		

For research use only.

GENE DESCRIPTION	chromosome 6 open reading frame 159
GENE ID	134701
ORIGIN SPECIES	Human
ACCESSIONS	BC130460
VECTOR	pDONR223
NCBI SYMBOL	C6orf159
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagacgcggaggcgcagaggtacagagagtgagctggcgtgcggccaccgacggccctacggcgcgcgggcgcggaactccgatacgaggtcttgagaccctgggtggacgccggaggaagaagaaggagagacgccgaaccagccgcgg
aggcgatcccgatggcctggaatgaccgagcctcaggaagctttaccaattcaggcaccagtcagactatftggccaaatcaaatgttatgattacttatcaagaagcagaagctctctgaaaaattttcaattcaagccacaatttcatttatgaagattctgtagcgaagat
gaaattgaggatctgacctgtgaaaattaa