

DEFB134-defensin, beta 134 Gene

Reference: PTXBC130461

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms:	beta-defensin 134; beta-defensin 135; defensin beta 134		

For research use only.

GENE DESCRIPTION	defensin, beta 134
GENE ID	613211
ORIGIN SPECIES	Human
ACCESSIONS	BC130461
VECTOR	pDONR223
NCBI SYMBOL	DEFB134
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagcctctcctgtgtgtttgtcttctttcccttgggatccagtctggcaggtataaattcattatcatcagaaatcacaaagaaatgctataaaaatggcatctgcagacttgaatgctatgagagtgaaatgttagttgcctactgtatgtttcagctggagtctgtgtcaaggaaatcctgcaccctga