

RP11-410N8.4-hypothetical LOC149950 Gene

Reference: PTXBC130462

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC149950		

For research use only.

GENE DESCRIPTION	hypothetical LOC149950
GENE ID	149950
ORIGIN SPECIES	Human
ACCESSIONS	BC130462
VECTOR	pDONR223
NCBI SYMBOL	RP11-410N8.4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggattgtaaatcccaaaagcgtgcaaacatctcccccatctccctggagggtggcctgtttctccacgccacctcacaggctgctgggggaactctctcacagcactgtgttcccggggcccactgttacaggccccatgcgagaaggacctcgggcagtgataatcctccaaggg
cccacaggaaacagcagtgacaactgtgttatgaagcactactatgtctgggggacaagaatggaacaagaagacatgctttgcctcaccctggagggcagcttcagccaggaggaataatctctccccagactcttccacagatccccaaactctgctctctcctcagatcc
ctgctcagatgtcaccagctctgtgcggcctcctgccaccaagtctaaaattgccatag