

LCN8-lipocalin 8 Gene

Reference: PTXBC130465

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EP17; LCN5; epididymal-specific lipocalin-8; lipocalin 5; lipocalin 8		

For research use only.

GENE DESCRIPTION	lipocalin 8
GENE ID	138307
ORIGIN SPECIES	Human
ACCESSIONS	BC130465
VECTOR	pDONR223
NCBI SYMBOL	LCN8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggagctggaccggcagaagattggaggattctggagggaagtcggtgtgacctcgatcaaacctgggtgtgacggccccgaagcgggtggaggcctgtgtctcacccttgagcgggagtaacctgacctgaaggttgcataaacagctcaggaagctgtgagatagaga
agatcgtgggtcagaaatagacagtacgggaaattcgttttctggtccacagagatccacgtgctggacaccgactacgagggtacgccatcctcgggtgtccctgatgtggcggggcaggaaactttcgctcctcaagtactttactcggagccttgaggacaaggaccggc
tggggttctggaagtttcgggagctgacagcagacactgggtctctacctggcggcccgccctgggcgggtgtgccgagctcctgaaggaggagctgatttaa