



Reference: PTXBC130468

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SPCA; coagulation factor VII; FVII coagulation protein; coagulation factor VII (serum prothrombin conversion accelerator); eptacog alfa; proconvertin		

For research use only.

GENE DESCRIPTION	coagulation factor VII (serum prothrombin conversion accelerator)
GENE ID	2155
ORIGIN SPECIES	Human
ACCESSIONS	BC130468
VECTOR	pDONR223
NCBI SYMBOL	F7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]