



Reference: PTXBC130471

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	F-box-like/WD repeat-containing protein TBL1Y; TBL1; transducin beta-like 1; transducin beta-like protein 1Y; transducin-beta-like protein 1, Y-linked; transducin beta like 1, Y-linked		

For research use only.

GENE DESCRIPTION	transducin (beta)-like 1Y-linked
GENE ID	90665
ORIGIN SPECIES	Human
ACCESSIONS	BC130471
VECTOR	pDONR223
NCBI SYMBOL	TBL1Y
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]