



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

KIF25-kinesin family member 25 Gene

Reference: PTXBC130489

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	kinesin-like protein KIF25; kinesin-like 3; kinesin-like protein 3; kinesin family member 25		

For research use only.

GENE DESCRIPTION	kinesin family member 25
GENE ID	3834
ORIGIN SPECIES	Human
ACCESSIONS	BC130489
VECTOR	pDONR223
NCBI SYMBOL	KIF25
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]