

IGFL2-IGF-like family member 2 Gene

Reference: PTXBC130490

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UNQ645; VPRI645; insulin growth factor-like family member 2; IGF like family member 2		

For research use only.

GENE DESCRIPTION	IGF-like family member 2
GENE ID	147920
ORIGIN SPECIES	Human
ACCESSIONS	BC130490
VECTOR	pDONR223
NCBI SYMBOL	IGFL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggaccgactacccagagtgctgctcctgcttatgtgtcagctgtctcctctgtgtccaagggaagtcacgctcccgctggctcagaacctggctgtccagccggcaccaggtgtggagacaagatctacaacccctggagcagtgctgttacaatgacgccatcgtgtccctgagcagagaccgccaatgtgtgtcccccgtgaccttctggccctgctttgagctctgctgtctttgattcctttggcctcacaacgattttgtgtgaagctgaaggttcagggtgtgaattccagtgccactatctccatctccagtaaatgtgaaagcagaagacgtttccctga