

C5orf48-chromosome 5 open reading frame 48 Gene

Reference: PTXBC130495

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C5orf48; Tseg7; testis-expressed sequence 43 protein; testis specific expressed gene 7; testis expressed 43		

For research use only.

GENE DESCRIPTION	chromosome 5 open reading frame 48
GENE ID	389320
ORIGIN SPECIES	Human
ACCESSIONS	BC130495
VECTOR	pDONR223
NCBI SYMBOL	C5orf48
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcttcagggaagatacttgcctacttgcctaaactactaacaactgctctgatgagagtcctataaatctgctaataagtatgaagattcattgcccagatttcattaagcaagggatgatccaagacgttatgtcatgccttgaaagaaaacatgatattcaggaatgtgaat
ctgaagcaagcagaagtgtgtggatccatactggcccttagaagactctctgttttgaatcacagtgaaggcctttgcccattggggaagatcgtaaagtgtcttccaaaaggcccaccagaaataaaattgcagatatgcctttgcattcgctctctccagatacacaagcactgtgattt
cccatggcttcaggaggcgactactaa