

KRTDAP-keratinocyte differentiation-associated protein Gene

Reference: PTXBC130501

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KDAP; UNQ467; keratinocyte differentiation-associated protein; KIPV467; keratinocyte differentiation associated protein		

For research use only.

GENE DESCRIPTION	keratinocyte differentiation-associated protein
GENE ID	388533
ORIGIN SPECIES	Human
ACCESSIONS	BC130501
VECTOR	pDONR223
NCBI SYMBOL	KRTDAP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagatccccgtctctgcccgtggtgctcctcctcctggtgctccactgcccaggagccaccctgggtggtcctgaggaagaaagcaccattgagaattatgcgtcacgacccgagcccttaacacccgttctgaacatcgacaattgcgactgcgtttaaggctgatg
agttcctgaactggcagccctctttgagctatcaaaaggaaactctcttctcaactgggatgcttcttaagctgaaggactgaggagcgaactcctgatgccagtgga