

**SPANXN4-SPANX family, member N4 Gene**

Reference: PTXBC130503

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT11.9; sperm protein associated with the nucleus on the X chromosome N4; cancer/testis antigen family 11, member 9; nuclear-associated protein SPAN-Xn4; SPANX family member N4		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	SPANX family, member N4
<b>GENE ID</b>	441525
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC130503
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	SPANXN4
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon )**

atggagagccaactccagcaccaacgagaataaaatgaagagccctgtgaatctaacaagaaaagttgacaagaagaagaatctgcacagagcctcagcccctgaacagagtttgaagagacagaaaaagcaaaatccaacattagttttactgcaggaagaataa  
gaaaaaaattcaaatcaactggagaataaccagcctacagagagctccactgatccaatcaagagaagagagacctagacatatctcaggatctccacaggatggtggcagaattag